


CUPS + SALIVA EJECTORS

Our cups' ribbed design offers a sure grip. Flexible saliva ejectors maintain their shape for effective hygiene.



TIDI® Choice 917902
Green 18" x 22" 3-ply tissue,
poly-back, contour neck bib.

9220 5-ounce TIDI® Tooth
print waxed paper cup.

PATIENT CUPS


Match colors and prints to our tissue and paper products for a friendly and coordinated exam environment.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9211	White	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9212	Green	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9213	Blue	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9214	Yellow	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9215	Gray	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9216	Mauve	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9217	Beige	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9218	Peach	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9210	Lavender	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9241	White	3.5 oz.	1000, 100/bag	Plastic cup, ribbed design
9243	Blue	3.5 oz.	1000, 100/bag	Plastic cup, ribbed design
9220	TIDI® Tooth	5 oz.	1000, 100/bag	Waxed paper cup
9221	TIDI® Tooth	5 oz.	500, 100/bag	Waxed paper cup
9226	Jazz	5 oz.	1000, 100/bag	Waxed infused paper cup
9240	Jazz	3 oz.	5000, 100/bag	Waxed infused paper cup
9231	White	6 oz.	2500, 250/bag	Evacuation cup liner, ADEC


COLORS



PATTERNS



Jazz



TIDI® Tooth



CUP DISPENSER

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9246	Chrome	3 or 5 oz., 16-1/4" x 3-3/4" x 3-3/4"	1	Adjustable cup dispenser

SALIVA EJECTORS

Saliva ejectors contain copper-plated wire, molded in, allowing them to bend easily and then hold their shape. Smooth tips are securely fastened.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9306	White	6"	1000, 100/bag	White ejector with white tip
9307	Clear	6"	1000, 100/bag	Clear ejector with blue tip
9304	White	Standard	1000, 100/bag	High-volume evacuator tip
9333	White	Standard	1000, 100/bag	High-volume evacuator tip, vented

